

theracool

COLD THERAPY SYSTEM

- Rheostat water flow control for temperature adjustment
- Integrated Pump System (IPS™)
- Quick-disconnect couplings
- Universal connector fits most pad brands
- Extremely quiet pump operation
- 9.6 quart capacity
- 90-day limited warranty (see policy)
- Specially designed circulation pads available for the hand, foot, shoulder, knee or in small and large sized multi-use configurations (sterile pads are available)
- Made in the USA



Theracool unit (CTS-1) and Slimline Knee Blanket (28591) sold separately

An Innovation in Cold Therapy

Sustained cold therapy is one of the most effective and soothing treatments for reducing post-operative pain and swelling. Well suited for hospital, home or portable use, the *TheraCool Cold Therapy System* helps return your patient to normal activity faster by providing cooling relief throughout the healing and recovery process. The TheraCool Home Unit allows patients to receive comforting, soothing cold therapy any time, anywhere.

Improving Circulation. The TheraCool exclusive pump lasts longer, and performs quietly. Its powerful pumping action helps to vary the level of cooling depending on patient need. An external integrated pump system helps reduce heat transfer. As a result, water stays cooler longer and leaves more space for ice and water.

The TheraCool Cold Therapy System. The TheraCool has a 9.6 quart capacity, easy-to-fill cooler with a secure lid to reduce accidental spills. A durable, long-lasting external pump and attached insulated hose assembly continuously circulate cold water through the specially designed, low-profile therapy blanket. Temperature is easily adjusted using a rheostat flow control which adjusts the flow of chilled water.

Your Cold Therapy Needs Are Covered. The *TheraCool Cold Therapy System* is available with a variety of Slimline Blankets designed to deliver continuous cold therapy to hip, knee, ankle, shoulder, back, or hand. The unique, unbroken loop covering on the outside of our Slimline Blankets allows for use without foam wraps. Simply use the provided elastic straps to secure the Slimline Blanket. They can even be used under post-operative bracing systems and are also available in sterile configurations (sterile blankets must be secured using foam wraps).

To get more information, or to locate a sales representative in your area, contact one of our customer service representatives today at (800) 222-4284.

Part#	Description	* CoolWrap required
CTS-1	Home Unit (universal connectors)	
28601	Slimline Blanket Hand/Foot	
28591	Slimline Blanket Knee	
28597	Slimline Blanket Shoulder	
28593	Slimline Blanket Small Rectangle [5"x10"]	
28595	Slimline Blanket Large Rectangle [8"x14"]	
28568-A	CoolPad Knee H ₂ O Circulation Pad*	
28567-A	CoolPad Shoulder H ₂ O Circulation Pad*	
28570-A	CoolPad Small Rectangle H ₂ O Circulation Pad* [5"x10"]	
28569-A	CoolPad Large Rectangle H ₂ O Circulation Pad* [8"x14"]	
28531-B	CoolWrap Knee Foam Wrap	
28551	CoolWrap Shoulder Foam Wrap	
28556	CoolWrap Small Rectangle Foam Wrap [5"x10"]	
28557	CoolWrap Large Rectangle Foam Wrap [8"x14"]	

Suggested HCPCS Codes: E0218 (unit), A9900 (pads)



visit us on the web at www.isports.com