

# GUARDIAN Sport EXT Rehabilitator™

Clinically more effective than  
*Dynamic Splinting and Static Progressive Stretch (SPS)*



- Patient preferred ambulating therapy reverses flexion contractures and extension deficits in gait
- Only therapy to simultaneously strengthen the quads while lengthening the hamstrings
- Patented extension assist technology provides an active extension stretch with every step
- Provides a therapeutic “glide” stretch with every step; estimated 90,000 active stretches in 30 days with 3 hours wear a day
- Eliminates Quadriceps Avoidance gait
- Ideal for post-surgery extension deficit therapy
- Ordered for knee instability due to extension deficit and generalized muscle weakness

BRACE SELECTION	INDICATION	STABILITY	PAIN RELIEF	UNBRACED PAIN RELIEF AFTER 90 DAYS	IMPROVE QUADS / HAMSTRING STRENGTH	IMPROVES KNEE EXTENSION	IMPROVES FOOT PLACEMENT	GAIT RETRAINING WITH RETAINED EFFECT	TREATS KNEE FLEXION CONTRACTURE
<b>SPORT EXT REHABILITATOR™</b>	Reverse Knee Flexion Contracture w/Gait	✓	✓	✓	✓	✓	✓	✓	✓
<i>Competitive Dynamic Brace</i>	Knee Flexion Contracture					✓			✓

Exclusively  
Provided by

**Bio DYNAMIC**  
TECHNOLOGIES

431 East 73rd Street • New York, NY 10021  
1 Madison Street, Suite A4 • East Rutherford, NJ 07073  
4 Dearfield Drive, Suite 201A • Greenwich, CT 06831

**To order call 1-800-879-2276  
or by fax 1-800-866-8011**